

## AtmoCheck<sup>®</sup> OPTIC O<sub>2</sub>

Optical Oxygen Analyser for random sampling of headspace



AtmoCheck<sup>®</sup> OPTIC O<sub>2</sub> is a laboratory analyser for the random quality control of residual oxygen concentrations. Utilising a state of the art optical oxygen sensor, our instrument provides reliable, fast and highly accurate results

The AtmoCheck<sup>®</sup> analyser works without a pump, and can therefore be used to accurately measure O<sub>2</sub> in extremely small headspace volumes. Examples of where this machine excels are: modified atmosphere packaging (MAP) with a minute headspace, or very small blister packs for pharmaceuticals.

AtmoCheck<sup>®</sup> analysers are in a class of their own and are characterised by decades of experience, exceptional measurement performance, robust design and simple intuitive operation.

In order to ensure and achieve the highest possible product quality, just rapidly test the residual oxygen content to the highest accuracy with AtmoCheck<sup>®</sup> OPTIC O<sub>2</sub>

### Food Industry Applications



The use of gases in the food industry in order to increase product quality and extend the shelf life of food is well established and helps to achieve:

- ➔ Quality Assurance
- ➔ Attractive Appearance
- ➔ Extended Shelf Life
- ➔ Fewer Preservatives

With our AtmoCheck<sup>®</sup> OPTIC O<sub>2</sub> analyser you can easily ensure the consistent high quality of your products by regularly checking the inert gas atmosphere. Ideal for quality control of MAP packaging with an extremely small headspace volume.

- |                    |                                    |
|--------------------|------------------------------------|
| ➔ MAP Packaging    | ➔ Coffee Packaging                 |
| ➔ Vacuum Packaging | ➔ Coffee Capsules, Coffee Pods     |
| ➔ Liquid Products  | ➔ Milk Powder, Milk Powder Sachets |

AtmoCheck<sup>®</sup> OPTIC O<sub>2</sub> is also optionally available for use directly on the packaging machine.

# AtmoCheck® OPTIC O<sub>2</sub>

## Pharmaceutical Industry Applications



Gases used in the production of active substances and pharmaceuticals are also used in a wide variety of areas and for various requirements, for example:

- ➔ as high purity gases for the laboratory
- ➔ as specialty gases for research and development
- ➔ as process gases for manufacturing processes such as chemical synthesis
- ➔ as inert gases for the removal of atmospheric oxygen
- ➔ as gases for sterile processes
- ➔ as traceable gases for GMP requirements
- ➔ as gas mixtures for measurements in the laboratory or in environmental
- ➔ as gases for the cultivation of organic crops

Many pharmaceutical products are very sensitive to atmospheric oxygen and moisture. Packaging in a protective atmosphere effectively prevents drugs from losing their effectiveness. The shelf life of therapeutic proteins in particular, which may be produced using nitrogen, is significantly extended.

AtmoCheck® OPTIC O<sub>2</sub> can optionally be used for the measurement of dissolved oxygen in a liquid.

## All highlights at a glance

- ➔ Application range for headspace volume < 2 ml
- ➔ Measuring range < 0.05% - 100%
- ➔ Short measuring time < 2 sec.
- ➔ Intuitive operation via 10" colour touchscreen display
- ➔ Integrated data logging software
- ➔ Graphic display
- ➔ Auto-Calibration Function
- ➔ Product menu and product management
- ➔ Ethernet connection for integration into networks
- ➔ Easy hygienic care due to splash-proof design
- ➔ Simple calibration in air via integrated SPAN procedure
- ➔ Low maintenance and robust design
- ➔ Also suitable for long-term testing

Your local partner will be happy to advise you and support you in the correct selection of analytical equipment. Our partners will also advise you on our individual or custom solutions for any requirement that you may have.

# AtmoCheck® OPTIC O<sub>2</sub>

## Special Features




AtmoCheck® OPTIC O<sub>2</sub> can be used in extreme conditions where traditional gas analysers cannot operate. The challenge of creating a sensor capable of measuring a minuscule headspace was to protect the sensor from food particles, powder or liquid. This was achieved with our new Optical Sensor Technology.

The outstanding feature of AtmoCheck® OPTIC O<sub>2</sub>, is that the Sensor is the only item that ever requires cleaning, and this is simplicity itself, taking literally only a few seconds.

In order to remove all soiling caused by either dry or liquid product residues that may remain on the sensor, all you need is just a glass of water. Simply place the tip of the sensor into the water and stir it for a couple of seconds. The sensor is then ready for use again within the shortest possible time.

The AtmoCheck® OPTIC O<sub>2</sub> is therefore virtually impervious to both solid and/or liquid contaminants.



AtmoCheck® OPTIC O<sub>2</sub> is intuitive, easy and safe to use. With a simple tap on the user interface you can control all functions; select products, create new products or users, etc.

Display of measured values, selected product, measuring time and limit values are clearly displayed on the large 10" touch display.

### Technical data

Gas	O <sub>2</sub> (residual oxygen) < 0.0 - 100%	Headspace volumes	< 2 ml
Measuring time	< 2 sek.	Interfaces	Ethernet, USB
Resolution	0,01% absolute *	Housing	Stainless steel
Accuracy	+/- 0,05% absolute **	Protection class	IP40
Service Life of O <sub>2</sub> Sensor	< 18 month	Weight	Approx. 7 kg
Heating time	< 1 minute	Dimensions	(HxWxD) 280 x 370 x 260 mm
Measuring temperature range	From 0° to + 50°C	Power supply	110 - 240V 50-60Hz. 50W

#### Manufacturing Accrediations and Standards

- Company certified according to ISO 9001:2015      - CE marked

\* Resolution:            ± 0,01 % O<sub>2</sub> at 1 % O<sub>2</sub>  
                                  ± 0,05 % O<sub>2</sub> at 20,9 % O<sub>2</sub>  
\*\* Accuracy at + 20 °C    ± 0,05 % O<sub>2</sub> or ± 3 % rel.